STEELMAN START

UNIVERSALUS NEIĞALIOJO VEŽIMĖLIS
UNIVERSĀLAIS RATIŅKRĒSLS
RATASTOOL
MANUAL WHEELCHAIR
ИНВАЛИДНАЯ КОЛЯСКА

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- Wheelchair is equipped with flip up and detachable armrest. To flip it up, please push plastic lever located in the front of armrest (2). Keep it pushed and lift armrest gently. We recomend to flip armrest up getting in or out of wheelchair from side. While operating wheelchair, armrests should be locked. In case the need to completely remove armrest occurs, firstly flip it up as explained before. Once armrest is flipped up, fnd plastic handle in the rear of armrest (3) and pull it out. Once pulled out, armrest is unlocked and you may take it out. To put it back, please repeate actions in apposite direction.

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
- To meet needs of different user, your wheelchair is equipped with height adjustable legrest. In order to adjust length of legrest, simply push metal button and pull the inned tube in or out depending you want to increase or decrease the length. You may set length in 6 different positions. The length of legrest should be set in position knee angle is bented in around 90°.

There is a risk of changing the point of gravity which leads to bad user's balance as well as risk of back bone deformation and bed sores development if you set length of legrest improperly.

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- Two different modes of rear wheel fixation are available in 04-020 series. The might be fixed with screw (version 04-020-1) or have quick release axle (1) (versions 04-020-2 and 04-020-3). If you have 04-020-1 version, you need to take a tool located in backrest pocket and loosen screw. Once loosened, hold frame with one hand and pull rear wheel out. To put it back, repeat action in apposite way. If your chair has quick release axle, you need no tool to take off rear wheel. Simply hold frame with one hand and push quick release axle located in the centre of rear wheel. Keep it pushed and pull rear wheel out. To put it back, repeat action in apposite way.

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OPERATION

For proper and safe use, we recommend you to follow safety instructions while using your wheelchair. This wheelchair was designed to ensure stability while seated. Every movement on the wheelchair affects the point of gravity, which may cause the wheelchair to tip over. To increase your safety, you may use seatbelt.

Reaching an object from the wheelchair - sitting

Only the arms may be extended over the wheelchairs seat. The body and head should stay inside the limits of the seat.

Lateral transfer

Transfer can be performed alone, assuming you are strong enough. Firstly, set the wheelchair as close to the place, which you want to sit on, as possible, whereby the steering wheels point forward. Activate the fixing brake and move your body towards the desired place. If there is no support beneath you while moving, please use a board if possible. Fold the legrest up and put them aside to make room. Bend your body forward, while leaning on the two armrests, and stand up. Then move your body toward the place, where you want to sit, while distributing your weight to your arms and hands.

Forward transfer

If you are capable of getting up and standing on your own, you may change seats by leaning forward. If you have trouble with this movement, we advise you to seek assistance from an accompanying person to ensure an easy and safe transfer. Firstly, engage both brakes and make sure both steering wheels point towards the front.

Please bear in mind that you must not support yourself on the brakes. Grab the handle of the armrest firmly and support yourself on your arms, to stand up easily (before, make sure armrest are securely locked in place). The assistant may help you complete the transfer. Position yourself as close to the desired place as possible. When changing your position, please sit as far back as possible, to prevent screws from breaking, seat cushions from getting damaged, or the wheelchair from tipping over. Activate the two brakes; these brakes should in no way be used as support for the transfer. You must never support yourself on the footrests, while sitting or getting up.

Moving onto the sidewalk

We recommend to have assistant to move onto sidewalk. Firstly, accompanying person positions the wheelchair straight in front of the curbside. After, tips the wheelchair backwards, until the steering wheels reach the curb. Lastly, pushes the wheelchair straight forward, until the back wheels hit the curb and when pushes the wheelchair again, until the back wheels go over the curbside.

Stairs

We suggest seeking help from two people, one in front of the wheelchair and one behind it, due to the fact that hurdling stairs is more difficult and dangerous. The person behind the wheelchair, tips the wheelchair backwards until it is balanced. Then, pushes the wheelchair against the first step, grabbing the handle firmly, to prepare to lift the wheelchair. The person in front of the wheelchair grabs the sides of the frame in the front of the wheelchair and lifts the chair over the steps, while the other accompanying person steps on the next step, repeating the first process. Do not lower the wheelchair until you hurdled the last step of the person in front of the wheelchair. Carrying a wheelchair down the stairs is the same process as explained previously, in reversed order.

Ascending a ramp

Lean your torso forward and use the hand rims with fast strokes, to keep up the proper speed. If you need a brake or relax on your journey, you must make sure to activate the two fixing brakes at the same time.

Descending a ramp

Lean your torso backwards and let the hand rims go (flow through your hands), making sure you have the speed and direction under control. Avoid sudden directional changes, never try to ascend or descend a ramp diagonally and never try to turn the wheelchair while ascending. The wheelchair is not applicable as a seat for cars.

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- We issue a warranty of 24 months from date of purchase. Any unauthorized modifications carried out without approval of the manufacturer will make the liability null and void. This includes exceeding the maximum user loads and using product in improper way. Warranty is only valid with fully filled warranty card. UAB „Kasko Group“ after getting claimed item from purchaser within period of 4 weeks decides whether issue is treated under warranty. If issue is not covered under warranty and customer agrees, UAB „Kasko Group“, if possible, fixes issue. All costs, including costs of transportation, should be covered by customer. Warranty repair does not extend warranty.

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